
ANDROMEDA VACUUM XP

MANUALE TECNICO / TECHNICAL MANUAL

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LIST OF CONTENTS

1. OPERATING CONTROLS AND DISPLAYS

- 1.1 FRONT VIEW – DOOR CLOSE
- 1.2 FRONT VIEW – DOOR OPEN
- 1.3 REAR VIEW

2. STERILIZATION CYCLES

- 2.1 OVERVIEW – CYCLES DESCRIPTION
 - 2.1.1 STERILIZATION CYCLES
 - 2.1.2 TEST CYCLES
- 2.2 STERILIZATION CYCLES TABLE
- 2.3 TEST CYCLES TABLE

3. UNPACKING AND INSTALLATION

4. SAFETY

- 4.1 SAFETY LABELS
- 4.2 SAFETY DEVICES
- 4.3 SAFETY NOTES

5. TECHNICAL DATA

6. ERROR MESSAGES

GUARANTEE


SPARE PARTS LIST

APPENDIXES (common for all languages)

- A. EXPLODED DRAWING
- B. FUNCTIONAL COMPONENTS
- C. PNEUMATIC DIAGRAM
- D. ELECTRIC DIAGRAMS
- E. CYCLE DIAGRAM
- F. PERIPHERALS TEST
- G. MAIN BOARD INFORMATION


1. OPERATING CONTROLS AND DISPLAYS

1.1. FRONT VIEW – DOOR CLOSE

- | | | | |
|----------|-----------------------|-----------|---------------------|
| 1 | Water level max | 8 | 134°C Special Cycle |
| 2 | Water level min | 9 | Cycle Select |
| 3 | Waste water level max | 10 | Water Inlet |
| 4 | 121°C Unwrapped Cycle | 11 | ON - OFF |
| 5 | 134°C Unwrapped Cycle | 12 | Open Door |
| 6 | 121°C Wrapped Cycle | 13 | Mains Switch |
| 7 | 134°C Wrapped Cycle | | |


1.2. FRONT VIEW – DOOR OPEN

- 1 Manual inlet water fitting
- 2 Microbiological Filter
- 3 Water inlet fitting
- 4 Waste water out

1.3. REAR VIEW

- 1 Discharge of used water
- 2 Exceeding used water – condensate discharge
- 3 Discharge clean water – connection deionizer
- 4 Exceeding clean water

- 5 Plug for Electrical power ONLY AND EXCLUSIVELY FOR DEIONIZER
- 6 Safety valve
- 7 Fuse 12A
- 8 Fuse 12A
- 9 Fuse 500mA
- 10 Connection RS232
- 11 Principal electrical alimentation

2. STERILIZATION CYCLES

2.1. OVERVIEW – CYCLES DESCRIPTION

The unit sterilizes and dries dental instruments according to the following cycles.

The instruments may be wrapped and then positioned on the trays into the sterilization chamber.

Sterilization cycles are classified in two categories, one of normally cycles used for sterilizing instruments, and the other for testing the autoclave.

2.1.1. STERILIZATION CYCLES

- **121°C unwrapped** is related to a solid load without wrapping
- **121°C wrapped** is related to a wrapped solid load
- **134°C unwrapped** is related to a solid load without wrapping
- **134°C wrapped** is related to a wrapped solid load
- **134°C Special** is related to a cycle with fast drying FLASH

2.2. STERILIZATION CYCLES TABLE

CYCLES	STERILIZ. TIME (Minutes)	DRYING TIME (Minutes)	MAX LOAD per Tray (Kg)	WORK PRESS. (bar)	VAC. PHASE N°
121°C unwrapped	18	17	0.75	1.04	1
121°C wrapped	21	17	0.75	1.04	1
134°C unwrapped	6	17	0.75	2.02	1
134°C wrapped	11	17	0.75	2.02	1
134°C special	6	2	0.75	2.02	1

3. UNPACKING AND INSTALLATION

Autoclave is supported with a wooden pallet and is inclosed in a corrugated board, inserted on the top of the machine.

The packaging is tied with two plastic straps to the pallet.

To unpackage the autoclave first cut the straps, than remove the board with a vertical movement.

Two persons are required to lift the machine, as its weight is 60 Kg.

Take care to the following notes:

- Install the autoclave only in rooms in which only authorized personnel is admitted
- The autoclaves do not have to be installed in places to risk upgrade them of outbreak and/or fire
- Autoclave room needs a good lightening
- Autoclave room needs a good ventilation
- Horizontal support plane load up to 70 Kg
- A 50 mm gap between the autoclave and the wall is needed
- Avoid the presence of water sprays towards the cover of the autoclave
- Do not put anything on the cover
- Do not install the machine near heat sources
- Avoid bending mains cable
- Take care that waste water tank level be at least 100 mm below machine bottom
- Take care that a.c. power supply be in compliance with safety rules.
- Take care that a.c. power supply be 230 V 2000 VA
- Do not insert mains connector through adapters

4. SAFETY

4.1. SAFETY LABELS

DANGEROUS VOLTAGE

DISCONNECT THE MAINS SUPPLY BEFORE REMOVING THIS COVER

HIGH TEMPERATURE

GROUND CONNECTION

4.2. SAFETY DEVICES

Safety devices are the following:

- safety pressure valve rated 2,4 Bar $\pm$ 10 %
- electromagnetic locking door to avoid opening the door when a cycle is running

4.3. SAFETY NOTES

As manufacturers of autoclaves we can only assume responsibility for the safety related performance of the unit, if **maintenance and repair** are carried out by us, or by agencies expressly authorized by us for this purpose, and if needed components affecting the safety are replaced exclusively by **original spare parts**.

The premises where it settles the autoclaves does not have to be to risk upgrades them of outbreak and/or fire.

In order to prevent injuries to persons and damage to the equipment you must take carefully into account all instruction that are emphasized with WARNING.

The autoclave must be installed at a distance at least 1,5 meters from the patient.

5. TECHNICAL DATA

5.1. MECHANIC

• Working temperature	+5 +40°C
• Max altitude	2.000 m
• Max relative humidity at 30°C	80%
• Max relative humidity at 40°C	50%
• Overall dimensions (wide x high x deep)	455x390x560 mm
• Open door space	370 mm
• Weight (empty tanks)	40 Kg
• Weight (full tanks)	43 Kg
• Pressure on bearing	1654 N\m ²
• Volume	0.12 m ³
• Colour	RAL 9016
• Materials	Stainless Steel AISI 304 Steel FeP01
• Noise	< 70 db A

5.2. ELECTRIC

• Supply Voltage	230 V a.c. + / - 10 % single phase
• Power	2.000 W
• Frequency	50 / 60 Hz
• Mains cable	2 + 1 x 1mm ²
• Fuses	6.3 x 32 - 12 A
• Heat production	3.6 E +6 J \ hour

5.3. CHAMBER

• Max working pressure	2.2 bar (relative)
• Max vacuum	- 0.9 bar (relative)
• Max temperature	138 °C
• Material	Stainless Steel AISI 304
• Dimensions	diam. 245 x 318 mm

5.4. CLEAN WATER TANK

• Volume	2.2 l
• Allowable cycles	5
• Material	Poliethylene

5.5. WASTE WATER TANK

- Volume 2.2 l
- Allowable cycles 5
- Material Polietylene
- Waste water max temp. 50°C

5.6. BACTERIOLOGICAL FILTER

- Diameter 56 mm
- Filtering efficiency 0.3 micrometer

5.7. TRAYS HOLDER

- Material Anodized Aluminium
- Dimensions (wide x high x deep) 192 x 165 x 280 mm
- Provided number 1

5.8. TRAYS

- Material Anodized Aluminium
- Dimensions (wide x high x deep) 184 x 17 x 286 mm
- Provided number 4

6. TECHNICAL DATA

Error	Time*	Problem	Solutions
E200	START	The cycle has been started with an insufficient level of water in the cargo tank	1. Refill the clean water tank 2. Check LEVEL1 3. Check connections between LEVEL1 and CBL607
E250	START	Level water in the cargo tank under the maximum level with the water network inlet	1. Wait for filling water tank 2. Check electrical and hydraulic connections of the water clean 3. Check LEVEL1 4. Check connections between LEVEL1 and CBL607
E300	START	The cycle has been started with the waste tank full	1. Empty the waste water tank 2. Check LEVEL2 3. Check connections between LEVEL2 and CBL607
E401	START	The cycle has been started with the door open	1. Chiudere il portello 2. Check SWITCH1 3. Check connections between SWITCH1 and CBL606B
E451	START	The cycle has been started with the electromagnet in the position "close"	1. Press the button "Open Door" per to unblock the safety pin 2. Check SWITCH1 and SWITCH2 3. Check the connections between SWITCH1, SWITCH2 and CBL606B
A001		Voluntary interruption of the sterilization cycle (caused by the pressure of the button I/O)	
A011		CBL606B doesn't communicate with CBL607	1. Check the connections between the two board and the state of both
A101	T4	Unable to complete the vacuum phase	1. Check the state of the door 2. Check EV2, EV3, EV4 3. Check the connections near the chamber 4. Check VALVE1 5. Check the brass filter
A501	T6	The chamber has lost the pressure	Like A101
A551		A temperature probe doesn't work	1. Check connections between TEMP1 and CBL606B 2. Replace TEMP1
A552		A temperature probe doesn't work	1. Check connections between TEMP2 and CBL606B 2. Replace TEMP2

A641	T5, T6	High pressure	1. Check PRESSURE1 2. Check connections between EV3, before T5
A651	T6	High temperature	1. Check EV2, EV4 2. Check the connections under the chamber 3. Check lose of water before T6 4. Check TEMP1, TEMP2
A661	T6	Temperature probes show different values	1. Check TEMP1, TEMP2
A701	T1, T3, T5	The pressure doesn't rise well	Like A101
A751	T6	Low temperature	1. Check TEMP1, TEMP2 2. Check HEAT1 3. Check lose of water or vapor
A761		Pressure sensor doesn't work correctly	1. Check connections between PRESSURE1 and CBL606B 2. Replace PRESSURE1
A801	T3, T4, T7	The pressure doesn't come down well	1. Check EV3, EV8 2. Check EV2
A901		Main supply voltage low	1. Check cable state

*** this mean the time interval in witch is possible to find the related error (see pressure diagram, APPENDIX E)**

GUARANTEE

1. The appliance is covered by a warranty period of twelve months.
2. The warranty begins the date the machine is delivered to the customer. This date must appear on the warranty coupon, which must be properly filled in, stamped and signed by the dealer and then returned to TECNO-GAZ.
3. In case of complaint, the date appearing on the bill of sale will be considered.
4. Repairs or replacement under warranty conditions of any components will be exclusively carried out by TECNO-GAZ; the warranty does not cover the travel expenses of personnel, packaging or transport expenses.
5. The warranty does not cover lamps, fuses, glass or failures or damage caused by incorrect maintenance, negligence, lack of expertise or any other cause the manufacturer may not be held not responsible for.
6. Those parts and equipment subject to normal wearing (e.g. buttons, pump moving parts, etc.) are not covered by the warranty.
7. TECNO-GAZ does not recognize any right the complete replacement of the machine.
8. The warranty does respond to any claims for damages, direct or indirect, to persons or things caused by the inefficiency of the machine.
9. TECNO-GAZ shall not be liable for failures or damage caused by incorrect use of the product or by failure to carry out ordinary maintenance or to observe the basic rules of its proper maintenance (negligence).
10. TECNO-GAZ does not recognize any right to indemnity in a case where the machine is not used.
11. The warranty will automatically fail to apply where the machine has been tampered with or repaired or modified by the client or third parties not authorized by TECNO-GAZ.
12. The purchaser must contact the dealer or after-sales services centre indicated by TECNO-GAZ for repairs, etc.
13. Components repaired during the warranty period must be returned, carriage paid, to TECNO-GAZ.
14. If replaced components are not returned to TECNO-GAZ the purchaser will be charged for them.
15. TECNO-GAZ does not admit restitution from final users.
16. Only the dealer or the after-sales assistance centre chosen by the user will be admitted to return to TECNO-GAZ components repaired according to commercial procedure CM-P-003 called "MANAGEMENT OF REQUESTS FOR THE RESTITUTION, REPAIR OR REPLACEMENT OF TECNO-GAZ EQUIPMENT AND CLIENT COMPLAINTS".
17. Components sent back to TECNO-GAZ must be accompanied by the correct documentation according to our internal procedures.
18. The products sent back to TECNO-GAZ must come with attached documentation authorising the return and a note describing the failure.
19. All components that need to be repaired must be sent to TECNO-GAZ carriage paid and properly packed (original packaging is recommended).

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

B. FUNCTIONAL COMPONENTS


Name	Description
CBL514D	Printer board
CBL606B	Power board
CBL607	Display board
EV2	EV normally close 1/4
EV3	EV normally close 1/4
EV4	EV normally close 1/8
EV5	EV normally open 1/4
EV8	EV normally close 1/8
FILTER1	Brass filter
FILTER2	Plastic filter
FILTER3	Bacteriological filter
PRESSURE1	Pressure sensor
VALVE1	Not return valve
VALVE2	Safety valve
MAINS	Mains switch
P230	Mains input 230V
FUSE1	Fuse 12A
FUSE2	Fuse 12A
PEXT	Mains auxiliary output (Water Clean)
FUSE3	Fuse 500mA
TANK1	Clean water tank
TANK2	Waste water tank
LEVEL1	Level clean water probe
LEVEL2	Level waste water probe
COOLER1	Cooler
FAN1	Fan
PUMP1	Load pump
PUMP2	Vacuum pump
VESSEL1	Vessel
BLOCK1	Door electromagnet
SWITCH1	Door open/close switch
SWITCH2	Door electromagnet switch
HEAT1	Chamber heat
TEMP1	Temperature probe 1
TEMP2	Temperature probe 1


TEMP3	Resetting manual thermostat – heater element
OUT-1	Waste water out
OUT-2	Condensed water out
OUT-3	Water tank out
OUT-4	Overflow out
OUT-A	Waste water out
IN-B	Water inlet fitting
IN-C	Manual inlet water fitting
RS232	RS232 connector
PANEL1	Front panel
PRINTER1	Printer

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

Pressure diagram when we press the key START for a cycle of 121°C or 134°C

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
F. PERIPHERALS TEST

It's possible to power up one by one all the peripherals connected to the main board without following the sterilization cycle. Attention: during the test, for example, if you power up the heating element then this will be always on until you don't stop it.

We suggest you to execute the test only with the door opened to avoid that the unit can go in pressure.

To run the test, proceed as follows. Insert **JUMPER** on **MAIN BOARD (DXGA058)** in **J3** position. Turn on the autoclave and move in setup menu' with **START** button to read up on the display USCITA. Use **PUMP WATER** button to increase the number on the display. At each number corresponds to the activation of a device connected to the card (see the table below) :

1. HEAT1 (!!! ATTENTION !!!)	-> resistance on
2. HEAT1 (!!! ATTENZIONE !!!)	-> resistance off

3. ---	
4. EV5	
5. EV2	
6. it closes relay "the direction" of BLOCK1	
7. BLOCK1	-> door closure 200VDC
8. ---	
9. it opens relay "the direction" of BLOCK1	
10. BLOCK1	-> opened door 200VDC
11. EV4	
12. PUMP1	-> water loading pump
13. PUMP2	-> vacuum pump
14. ---	
15. EV3	
16. ---	
17. ---	
18. EV8	
19. BUZZER	

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CN14 = Collegamento a scheda display

G. MAIN BOARD INFORMATION

The DXGA058, is the main board of the autoclave. It rules :

-) heating band
-) solenoid
-) vacuum pump
-) water loading pump
-) door closer system

IMPORTANT

Before installing the board verify the correspondence of software. In the case, Remove the microprocessor from the old board and insert the new one.

Attention to this operation because very delicate.

Once inserted the microprocessor must enter the serial number of the machine and the number of the customer.


Numbers to be inserted
relating to the serial
number of the machine

Label example

IMPORTANT

annotate also inside the machine on the reservoirs the date of this intervention and the number of cycles of the machine.

CONNECTIONS DESCRIPTION

EPROM = Identifies the software, the dot is pin number 1

Fuses F1 e F2 = Protect all the services to the exclusion of the resistance

BATTERY = is the buffer battery

J3 = Jumper to insert to enter in the technical test

CN1 = Collegamento Servizi

Pin 1 = Ground

Pin 2 = Charging pump water

Pin 3 = Vacuum pump

Pin 4 = ---

Pin 5 = EV 3

Pin 6 = Socket osmosi s system

Pin 7 = EV 4

Pin 8 = ---

CN2 = Services connection

Pin 1 = Ground (connected to PIN 1 with CN1)

Pin 2 = EV8

Pin 3 = EV2

Pin 4 = EV5

Pin 5 =

CN3 = Block1 electromagnet opening-closing door (is fed to 200 Vdc)

CN4 = Printer

CN9 – CN10 = Connection of the two probes that are located inside the chamber, the board is calibrated so that the difference of reading between the two probes is 0.01 °C, in case of alarms inherent to the malfunctioning of the probes try to reverse the probes connection between them.

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MEDILINE ITALIA S.r.l.
TECNO-GAZ S.p.A.